



# SAFE life®


## Fentanyl (FTY) Liquid & Powder Drug Test Strip

- Fentanyl Harm Reduction
- Rapid Results
- Accurate

SAFElife™ Fentanyl (FTY) Liquid & Powder Drug Test Strip is an immunochromatographic assay for the qualitative determination of Fentanyl in liquid, powder or pill samples.

### Background

Fentanyl (FTY) is a potent, synthetic opioid analgesic with a rapid onset and short duration of action. Fentanyl is approximately 50 times stronger than heroin and 100 times stronger than morphine. It has become a major contributor to fatal and nonfatal overdoses in the U.S. While it is used clinically for severe pain, most fentanyl-related overdose cases are linked to its presence in illegal drug markets. In response to this crisis, the **SAFElife™ Fentanyl (FTY) Liquid & Powder Drug Test Strip** can determine if fentanyl is present in a liquid or powder, providing vital information to help with harm reduction.



## Features and Benefits

- Detects fentanyl and multiple fentanyl analogs
- No false positives with MDMA, Methamphetamine, Cocaine or Morphine
- Results in 5 minutes
- 99% Accurate

Drug (Identifier)	Fentanyl (FTY)
Calibrator	Fentanyl
Cut-off Level	10 ng/mL
Sample Types	Liquid or Powder

## Test Procedure

**Step 1 Add water**

**Step 2 Test**

**Step 3 Wait**

**Step 4 Results**

Item Code	Description	Box Size / Case Size	Complexity
WSPFTY-101	1-Strip; Drug Liquid / Powder Detection; Fentanyl 10 ng/mL	25/120	Forensic Use Only Harm Reduction Use

Rel.: 2023/11/13