




SAFE life™

T-Swab Multi-Drug Oral Fluid Test



- Tests can be run within 2 minutes
- Easy to use
- E&I Exempt
- Screens for 7 common drugs of misuse
- Result at 5 minutes

Employment and Insurance Testing:

The Food and Drug Administration (FDA or Agency) is announcing a list of class II devices that the Agency has determined based on established factors to no longer require premarket notification to provide reasonable assurance of safety and effectiveness, subject to certain limitations. FDA is publishing this notice of that determination in accordance with procedures established by the 21st Century Cures Act.

Workplace Drug Testing:

Follow-Up, for cause, periodic, post-accident, pre-employment, random, return to duty. Not applicable for Federally regulated workplace and DOT testing programs.

T-Swab™ Cutoff Levels

Drug Name	Cut off (ng/ml)
AMP (Amphetamine)	50
COC (Cocaine)	20
THC (Marijuana)	40
MET (Methamphetamine)	50
OPI (Opiate)	40
OXY (Oxycodone)	20
PCP (Phencyclidine)	10

We offers 4 combinations from 7 substances.

Item Code	Test Combination
** SODOA-156NTEUO	AMP-COC-OPI; MET-PCP
** SODOA-256EUO	AMP-COC-OPI; MET-THC
SODOA-166EUO	AMP-COC-OPI; MET-THC-PCP
** SODOA-666EUO	AMP-COC-OPI; MET-THC-OXY

** Ordered upon request