



WELLlife™

COVID-19 / Influenza A&B Test

The WELLlife™ COVID-19 / Influenza A&B Test differentiates between SARS-COV-2, influenza A and influenza B antigens with a single test.

Features and Benefits

- **Quality & Accuracy**

510(k) cleared for professional use.

- **Cost-effective**

Differentiate between SARS-COV-2, influenza A, and influenza B antigens with a single test using an anterior nasal swab specimen.

- **Rapid**

Results available in 10 minutes allow for testing and treatment decision-making during the same office visit.

- **Quality Assurance**

The outstanding clinical performance, including built-in kit controls for external quality testing, enables physicians to make quicker and more confident decisions.

- **Extended Detection Window**

Offers a broader detection window for differentiation of SARS-COV-2, influenza A, and influenza B antigens within four (4) days of symptom onset.

- **Wide Storage Temperature**

36°F–86°F (2°C–30°C) allows for easier product storage.



Ordering Information

Item Code	Description	Size	Intended Use
WV01P0002	WELLlife™ COVID-19 / Influenza A&B Test	25T/Kit	For professional use

Rel.: 2025 / 01 / 31