



SAFE life™

Fentanyl (FTY) Liquid & Powder Drug Test Strip

- Fentanyl Harm Reduction
- Rapid Results
- Accurate

SAFElife™ Fentanyl (FTY) Liquid & Powder Drug Test Strip is an immunochromatographic assay for the qualitative determination of Fentanyl in urine, liquid, powder or pill samples.

Background

Fentanyl (FTY) is a potent, synthetic opioid analgesic with a rapid onset and short duration of action. Fentanyl is approximately 50 times stronger than heroin and 100 times stronger than morphine. It has become a major contributor to fatal and nonfatal overdoses in the U.S. While it is used clinically for severe pain, most fentanyl-related overdose cases are linked to its presence in illegal drug markets. In response to this crisis, the **SAFElife™ Fentanyl (FTY) Liquid & Powder Drug Test Strip** can determine if fentanyl is present in a urine, liquid, powder or pill, providing vital information to help with harm reduction.




Features and Benefits


- Detects fentanyl and multiple fentanyl analogs.
- No false positives with MDMA, Methamphetamine, Cocaine or Morphine.
- Results in 1 minute.
- Easy to Use: Each strip features clear MAX Line and MIN Line indicators.
- Cost-effective option: Test Strip Version and Accessory Kit Version to meet diverse needs.
- All-in-One Convenience: All the accessories you need in one kit for user-friendly simplicity.

Drug (Identifier)	Fentanyl (FTY)
Calibrator	Fentanyl
Cut-off Level	10 ng/mL
Reading Time	1 min
Sample Type	urine, liquid, powder or pill


Detection of Fentanyl and Its Analogues			
Carentanil	Norfentanyl	Isobutyryl fentanyl	3-Methylfentanyl
Acryl fentanyl	p-Fluorofentanyl	4-Fluoro-isobutyryl fentanyl	Valeryl fentanyl
Cyclopropyl fentanyl	Butyryl fentanyl	para-Fluorobutyryl fentanyl (PFBF)	Ocfentanil
Methoxyacetyl fentanyl	Furanyl fentanyl	Acetyl fentanyl	o-Fluorofentanyl
Remifentanil			



Liquid & Powder Drug Test Strip



Liquid & Powder Drug Test Kit Pro




Liquid & Powder Drug Test Kit Premier




SAFE life™

Test Procedure

Step 1 Add water




Step 2 Test



Step 3 Wait



Step 4 Result



Item Code	Description	Packing	Complexity
WSPFTY-101	1-Strip Drug Liquid / Powder Detection Test; FTY Fentanyl 10 ng/ml	25T	Forensic Use Only Harm Reduction Use
WSPFTY-101-IPAP	1-Strip Drug Liquid / Powder Detection with Mixing Accessories Pro Version; FTY Fentanyl 10 ng/ml	1 Test per Pack	Forensic Use Only Harm Reduction Use