



SAFE life™

T-Cup Multi-Drug Urine Test Cup

Rapid Point-of-Care Toxicology Testing

- CLIA Waived, 510(K) Cleared & Forensic Use Only
- Female User Friendly & Customizable
- Test up to 26 Drugs per Cup



CLIA-waived Drugs & Cutoff Levels

AMP AMP AMP	Amphetamine (AMP)	1000, 500 ng/mL
BAR BAR BAR	Barbiturates (BAR)	300 ng/mL
BZO BZO BZO	Benzodiazepines (BZO)	300 ng/mL
BUP BUP BUP	Buprenorphine (BUP)	10 ng/mL
COC COC COC	Cocaine (COC)	300, 150 ng/mL
THC THC THC	Cannabinoids (THC)	50 ng/mL
EDDP EDDP EDDP	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300 ng/mL
MTD MTD MTD	Methadone (MTD)	300 ng/mL

TCA TCA TCA	Tricyclic Antidepressants (TCA)	1000 ng/mL
MET MET MET	Methamphetamine (MET/mAMP)	1000, 500 ng/mL
MDMA MDMA MDMA	Methylenedioxymethamphetamine (MDMA)	500 ng/mL
MOP MOP MOP	Morphine (MOP/OPI300)	300 ng/mL
OPI OPI OPI	Opiate (OPI)	2000 ng/mL
OXY OXY OXY	Oxycodone (OXY)	100 ng/mL
PCP PCP PCP	Phencyclidine (PCP)	25 ng/mL
PPX PPX PPX	Propoxyphene (PPX)	300 ng/mL

Forensic Use Panels

6-MAM 6-MAM 6-MAM	Monoacetylmorphine (6-MAM)	10 ng/mL
AMP AMP AMP	Amphetamine (AMP)	300, 500 ng/mL
BUP BUP BUP	Buprenorphine (BUP)	5 ng/mL
BZO BZO BZO	Benzodiazepines (BZO)	100, 200 ng/mL
COC COC COC	Cocaine (COC)	100, 150 ng/mL
COT COT COT	Cotinine (COT)	200 ng/mL
EDDP EDDP EDDP	EDDP	100 ng/mL
EtG EtG EtG	Ethyl Glucuronide (EtG)	500 ng/mL
FTY FTY FTY	Fentanyl (FTY)	20 ng/mL
GAB GAB GAB	Gabapentin (GAB)	2000 ng/mL
HMO HMO HMO	Hydromorphone (HMO)	300 ng/mL

K2 K2 K2	Synthetic Cannabinoid (K2)	50 ng/mL
KET KET KET	Ketamine (KET)	300, 1000 ng/mL
KRA KRA KRA	Kratom (KRA)	100, 300 ng/mL
LSD LSD LSD	Lysergic Acid Diethylamide (LSD)	20 ng/mL
MET MET MET	Methamphetamine (MET/mAMP)	300, 500 ng/mL
MOP MOP MOP	Morphine (MOP or OPI)	100 ng/mL
MTD MTD MTD	Methadone (MTD)	200 ng/mL
MQL MQL MQL	Methaqualone (MQL)	300 ng/mL
PGB PGB PGB	Pregabalin (PGB)	500 ng/mL
THC THC THC	Cannabinoids (THC)	15, 25, 40 ng/mL
TRA TRA TRA	Tramadol (TRA)	100, 200, 1000 ng/mL

Generic stock adulteration testings include Creatinine (CR), Specific Gravity (S.G.) and pH. Other adulterates like Glutaraldehyde (GL), Oxidant (OX) & Nitrite (NI) are available through special order.



SAFE life™

Product Information

※ Kit Size: 25 / Case Size: 4

Item Code	Description	Complexity
TDOA-154	5-Panel T-Cup; AMP, COC, OPI, PCP, THC	CLIA WAIVED
TDOA-254	5-Panel T-Cup; AMP, COC, mAMP, OPI, THC	CLIA WAIVED
TDOA-254A3	5-Panel T-Cup; AMP, COC, mAMP, OPI, THC with 3 AD	CLIA WAIVED
TDOA-164	6-Panel T-Cup; AMP, COC, MET, OPI, PCP, THC	CLIA WAIVED
TDOA-264	6-Panel T-Cup; AMP, BZO, COC, mAMP, OPI, THC	CLIA WAIVED
TDOA-274	7-Panel T-Cup; AMP, BZO, COC, mAMP, OPI, OXY, THC	CLIA WAIVED
TDOA-3104	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, PCP, THC	CLIA WAIVED
TDOA-3104A3	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, PCP, THC with 3 AD	CLIA WAIVED
TDOA-4104	10-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MTD, OPI, OXY, PCP, THC	CLIA WAIVED
TDOA-7104	10-Panel T-Cup; BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, TCA, THC	CLIA WAIVED
TDOA-8104	10-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, THC	CLIA WAIVED
TDOA-4124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, PCP, PPX, THC	CLIA WAIVED
TDOA-6124	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC	CLIA WAIVED
TDOA-6124A3	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC with 3 AD	CLIA WAIVED
TDOA-6125	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA WAIVED
TDOA-6125A3	12-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC with 3 AD	CLIA WAIVED
TDOA-6125V	12-Panel T-Cup w/Vacu Lid; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,THC	CLIA WAIVED
TDOA-7125	12-Panel T-Cup; AMP, BAR, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC	CLIA WAIVED
TDOA-2135	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC	CLIA WAIVED
TDOA-2135A3	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC with 3 AD	CLIA WAIVED
TDOA-6135F	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, FTY	Forensic Use Only
TDOA-7135EFT	13-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG, FTY, TRA	Forensic Use Only
TDOA-7135A3EFT	13-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MTD, MOP, OXY, THC, ETG, FTY, TRA with 3 AD	Forensic Use Only
TDOA-8135E	13-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG	Forensic Use Only
TDOA-1144	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC	CLIA WAIVED
TDOA-1144A3	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, TCA, THC with 3 AD	CLIA WAIVED
TDOA-1145A3	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, PPX, TCA, THC with 3 AD	CLIA WAIVED
TDOA-8145EF	14-Panel T-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FTY	Forensic Use Only
TDOA-9145EFTK	14-Panel T-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG, FTY, TRA, K2	Forensic Use Only
TDOAC-L1223A5EFTKK6	22-Panel T-Cup Clear; AMP300,BAR300,BUP10,BZO200,COC100,EDDP300,HMO300,KET1000,MDMA500,MET300,MOP100,MTD200,OXY100,PCP25,PPX300,THC25,ETG500,FTY20,TRA200,K2 50,KRA300,6-MAM10 +(CR,NI,OX,pH,S.G.)	Forensic Use Only
FDOA-6125	12-Panel F-Cup; AMP, BAR,BUP,BZO,COC,mAMP, MDMA, MOP, MTD, OXY, PCP,THC	CLIA WAIVED

Rel.: 2024/05/17