




SAFE life™

# T-Dip Multi-Drug Urine Test Panel

## Rapid Point-of-Care Toxicology Testing

- CLIA Waived, 510(K) Cleared & Forensic Use Only
- Cost Effective & Customizable
- Test up to 16 Drugs per Card



## CLIA-waived Drugs & Cutoff Levels

|                |                                                          |                 |                |                                      |                 |
|----------------|----------------------------------------------------------|-----------------|----------------|--------------------------------------|-----------------|
| AMP AMP AMP    | Amphetamine (AMP)                                        | 1000, 500 ng/mL | TCA TCA TCA    | Tricyclic Antidepressants (TCA)      | 1000 ng/mL      |
| BAR BAR BAR    | Barbiturates (BAR)                                       | 300 ng/mL       | MET MET MET    | Methamphetamine (MET/mAMP)           | 1000, 500 ng/mL |
| BZO BZO BZO    | Benzodiazepines (BZO)                                    | 300 ng/mL       | MDMA MDMA MDMA | Methylenedioxymethamphetamine (MDMA) | 500 ng/mL       |
| BUP BUP BUP    | Buprenorphine (BUP)                                      | 10 ng/mL        | MOP MOP MOP    | Morphine (MOP/OPI300)                | 300 ng/mL       |
| COC COC COC    | Cocaine (COC)                                            | 300, 150 ng/mL  | OPI OPI OPI    | Opiate (OPI)                         | 2000 ng/mL      |
| THC THC THC    | Cannabinoids (THC)                                       | 50 ng/mL        | OXY OXY OXY    | Oxycodone (OXY)                      | 100 ng/mL       |
| EDDP EDDP EDDP | 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP) | 300 ng/mL       | PCP PCP PCP    | Phencyclidine (PCP)                  | 25 ng/mL        |
| MTD MTD MTD    | Methadone (MTD)                                          | 300 ng/mL       | PPX PPX PPX    | Propoxyphene (PPX)                   | 300 ng/mL       |

## Forensic Use Panels

|                   |                            |                |             |                                  |                      |
|-------------------|----------------------------|----------------|-------------|----------------------------------|----------------------|
| 6-MAM 6-MAM 6-MAM | Monoacetylmorphine (6-MAM) | 10 ng/mL       | K2 K2 K2    | Synthetic Cannabinoid (K2)       | 50 ng/mL             |
| AMP AMP AMP       | Amphetamine (AMP)          | 300, 500 ng/mL | KET KET KET | Ketamine (KET)                   | 300, 1000 ng/mL      |
| BUP BUP BUP       | Buprenorphine (BUP)        | 5 ng/mL        | KRA KRA KRA | Kratom (KRA)                     | 100, 300 ng/mL       |
| BZO BZO BZO       | Benzodiazepines (BZO)      | 100, 200 ng/mL | LSD LSD LSD | Lysergic Acid Diethylamide (LSD) | 20 ng/mL             |
| COC COC COC       | Cocaine (COC)              | 100, 150 ng/mL | MET MET MET | Methamphetamine (MET/mAMP)       | 300, 500 ng/mL       |
| COT COT COT       | Cotinine (COT)             | 200 ng/mL      | MOP MOP MOP | Morphine (MOP or OPI)            | 100 ng/mL            |
| EDDP EDDP EDDP    | EDDP                       | 100 ng/mL      | MTD MTD MTD | Methadone (MTD)                  | 200 ng/mL            |
| EtG EtG EtG       | Ethyl Glucuronide (EtG)    | 500 ng/mL      | MQL MQL MQL | Methaqualone (MQL)               | 300 ng/mL            |
| FTY FTY FTY       | Fentanyl (FTY)             | 20 ng/mL       | PGB PGB PGB | Pregabalin (PGB)                 | 500 ng/mL            |
| GAB GAB GAB       | Gabapentin (GAB)           | 2000 ng/mL     | THC THC THC | Cannabinoids (THC)               | 15, 25, 40 ng/mL     |
| HMO HMO HMO       | Hydromorphone (HMO)        | 300 ng/mL      | TRA TRA TRA | Tramadol (TRA)                   | 100, 200, 1000 ng/mL |

Generic stock adulteration testings include Creatinine (CR), Specific Gravity (S.G.) and pH. Other adulterates like Glutaraldehyde (GL), Oxidant (OX) & Nitrite (NI) are available through special order.



# SAFE life™

## T-Dip Single Drug Urine Test Card

※ Kit Size: 25 / Case Size: 20

| Item Code   | Description                                   | Complexity        |
|-------------|-----------------------------------------------|-------------------|
| WDBU-114    | 1-Panel Card; BUP Buprenorphine               | CLIA WAIVED       |
| WDBZ-114    | 1-Panel Card; BZO Benzodiazepines             | CLIA WAIVED       |
| WDCO-114    | 1-Panel Card; COC Cocaine                     | CLIA WAIVED       |
| WDMA-114    | 1-Panel Card; mAMP/MET Methamphetamine        | CLIA WAIVED       |
| WDMT-114    | 1-Panel Card; MTD Methadone                   | CLIA WAIVED       |
| WDOP-114    | 1-Panel Card; OPI Opiate 2000 ng/mL           | CLIA WAIVED       |
| WDOX-114    | 1-Panel Card; OXY Oxycodone                   | CLIA WAIVED       |
| WDPC-114    | 1-Panel Card; PCP Phencyclidine               | CLIA WAIVED       |
| WDTH-114    | 1-Panel Card; THC Marijuana/Cannabinoids      | CLIA WAIVED       |
| W6AM-114    | 1-Panel Card; 6-MAM Monoacetylmorphine        | Forensic Use Only |
| WCOT-114    | 1-Panel Card; COT Cotinine                    | Forensic Use Only |
| WDK2-114    | 1-Panel Card; K2 Synthetic Cannabinoid        | Forensic Use Only |
| ** WDTH-15  | 1-Panel Card; THC Marijuana 15 ng/mL          | Forensic Use Only |
| WDTR-114    | 1-Panel Card; TRA Tramadol                    | Forensic Use Only |
| WETG-114    | 1-Panel Card; EtG Ethyl Glucuronide           | Forensic Use Only |
| WETG-300    | 1-Panel Card; EtG Ethyl Glucuronide 300 ng/mL | Forensic Use Only |
| WFTY-114    | 1-Panel Card; FTY Fentanyl                    | Forensic Use Only |
| ** WGAB-114 | 1-Panel Card; GAB Gabapentin                  | Forensic Use Only |
| ** WHMO-114 | 1-Panel Card; HMO Hydromorphone               | Forensic Use Only |
| WKET-114    | 1-Panel Card; KET Ketamine 1000 ng/mL         | Forensic Use Only |
| WKRA-114    | 1-Panel Card; KRA Kratom                      | Forensic Use Only |
| WMQL-114    | 1-Panel Card; MQL Methaqualone                | Forensic Use Only |
| ** WPGB-114 | 1-Panel Card; PGB Pregabalin                  | Forensic Use Only |

\*\* Ordered upon request

## T-Dip Multi-Drug Urine Test Panel

※ Kit Size: 25 / Case Size: 20

| Item Code   | Description                                                                           | Complexity        |
|-------------|---------------------------------------------------------------------------------------|-------------------|
| WDOA-234    | 3-Panel Card; COC,mAMP,THC                                                            | CLIA WAIVED       |
| WDOA-154    | 5-Panel Card; AMP,COC,OPI,PCP,THC                                                     | CLIA WAIVED       |
| WDOA-254    | 5-Panel Card; AMP,COC,mAMP,OPI,THC                                                    | CLIA WAIVED       |
| WDOA-654    | 5-Panel Card; BZO,COC,mAMP,OPI,THC                                                    | CLIA WAIVED       |
| WDOA-754    | 5-Panel Card; AMP,BZO,COC,OPI,THC                                                     | CLIA WAIVED       |
| WDOA-264    | 6-Panel Card; AMP,BZO,COC,mAMP,OPI,THC                                                | CLIA WAIVED       |
| WDOA-564    | 6-Panel Card; BZO,COC,mAMP,OPI,OXY,THC                                                | CLIA WAIVED       |
| WDOA-865NT  | 6-Panel Card without THC; BUP,BZO,COC,MOP,MTD,OXY                                     | CLIA WAIVED       |
| WDOA-274    | 7-Panel Card; AMP,BZO,COC,mAMP,OPI,OXY,THC                                            | CLIA WAIVED       |
| WDOA-295    | 9-Panel Card; AMP,BUP,BZO,COC,MOP,MTD,OXY,PPX,THC                                     | CLIA WAIVED       |
| WDOA-994NT  | 9-Panel Card without THC; AMP,BAR,BUP,BZO,COC,MET,MTD,OPI,OXY                         | CLIA WAIVED       |
| WDOA-1104   | 10-Panel Card; AMP,BAR,BZO,COC,mAMP,MTD,OPI,PCP,TCA,THC                               | CLIA WAIVED       |
| WDOA-3104   | 10-Panel Card; AMP,BAR,BZO,COC,mAMP,MDMA,MTD,OPI,PCP,THC                              | CLIA WAIVED       |
| WDOA-4104   | 10-Panel Card; AMP,BAR,BZO,COC,mAMP,MTD,OPI,OXY,PCP,THC                               | CLIA WAIVED       |
| WDOA-7104   | 10-Panel Card; BAR,BUP,BZO,COC,mAMP,MTD,OPI,OXY,TCA,THC                               | CLIA WAIVED       |
| WDOA-9104   | 10-Panel Card; AMP,BAR,BUP,BZO,COC,mAMP,MTD,OPI,OXY,THC                               | CLIA WAIVED       |
| WDOA-1124   | 12-Panel Card; AMP,BAR,BZO,COC,mAMP,MDMA,MTD,OPI,OXY,PCP,PPX,THC                      | CLIA WAIVED       |
| WDOA-3124   | 12-Panel Card; AMP,BAR,BUP,BZO,COC,MDMA,MTD,TCA,OPI,OXY,PCP,THC                       | CLIA WAIVED       |
| WDOA-6124   | 12-Panel Card; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MTD,OPI,OXY,PCP,THC                      | CLIA WAIVED       |
| WDOA-6125   | 12-Panel Card; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,THC                      | CLIA WAIVED       |
| WDOA-7125   | 12-Panel Card; AMP,BAR,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,TCA,THC                      | CLIA WAIVED       |
| WDOA-8145EF | 14-Panel Card; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FTY | Forensic Use Only |

Rel.: 2024/05/10