



SAFE life™

T-Square Multi-Drug Oral Fluid Test

Rapid Point-of-Care Toxicology Testing

SAFElife's T-Square Multi-Drug Oral Fluid Test is a simple and convenient method for drug screening. Full observation during collection allows testing to be performed anytime and anywhere, which greatly reduces specimen tampering concerns. Because it takes time for drugs to metabolize and pass through the body, oral fluid is also recognized as the optimal method for impairment detection.

Key Benefits:

- Convenient rapid screening
- Test up to 12 drugs & alcohol
- E&I exempt and Forensic Use Only testing panels
- Rapid results
- Easily witnessed collection
- Tamper resistant
- Detects recent drug use



Expanded Application of Oral Fluid Drug Testing to Non-Federal Workplace and Insurance

The Food and Drug Administration (FDA or Agency) is announcing a list of class II devices that the Agency has determined based on established factors to no longer require premarket notification to provide reasonable assurance of safety and effectiveness, subject to certain limitations. FDA is publishing this notice of that determination in accordance with procedures established by the 21st Century Cures Act.

*The Employment & Insurance (E&I) Exemption Notice by the Food and Drug Administration on 07/11/2017.



SAFE life™

E&I Exempt Panels

Drug	Cutoff	Drug	Cutoff
Amphetamine (AMP)	50 ng/mL	Methamphetamine (MET)	50 ng/mL
Barbiturates (BAR)	20 ng/mL	Methylenedioxy-methamphetamine (MDMA)	50 ng/mL
Cocaine (COC)	20 ng/mL	Opiate (OPI)	40 ng/mL
Marijuana (THC)	40 ng/mL	Oxycodone (OXY)	20 ng/mL
Methadone (MTD)	30 ng/mL	Phencyclidine (PCP)	10 ng/mL

Forensic Use Panels

 Alcohol (ACL)	0.02 g/dL	 Ketamine (KET)	100 ng/mL
 Amphetamine (AMP)	50 ng/mL	 Methadone (MTD)	30 ng/mL
 Barbiturates (BAR)	60 ng/mL	 Methamphetamine (MET/mAMP)	50 ng/mL
 Benzodiazepines (BZO)	30 ng/mL	 Methylenedioxymethamphetamine (MDMA)	100 ng/mL
 Buprenorphine (BUP)	5 ng/mL	 Opiate (OPI)	40 ng/mL
 Cannabinoids (THC)	25 ng/mL	 Oxycodone (OXY)	20 ng/mL
 Cocaine (COC)	20 ng/mL	 Phencyclidine (PCP)	10 ng/mL
 Fentanyl (FTY)	30 ng/mL	 Marijuana (THC)	15 ng/mL

Product Information

※ Kit Size: 25 / Case Size: 20

Item Code	Description	Complexity
QODOA-156NTEUO-I	5-Panel T-Square no THC with indicator; AMP, COC, MET, OPI, PCP	E&I Exempt
QODOA-256EUO-I	5-Panel T-Square with indicator; AMP, COC, MET, OPI, THC	E&I Exempt
QODOA-166EUO-I	6-Panel T-Square with indicator; AMP, COC, MET, OPI, PCP, THC	E&I Exempt
QODOA-366NTEUO-I	6-Panel T-Square no THC with Indicator; AMP, COC, MET, OPI, OXY, PCP	E&I Exempt
QODOA-666EUO-I	6-Panel T-Square with indicator; AMP, COC, MET, OPI, OXY, THC	E&I Exempt
QODOA-376EUO-I	7-Panel T-Square with Indicator; AMP, COC, MET, OPI, OXY, PCP, THC	E&I Exempt
QODOA-496NTEUO-I	9-Panel T-Square no THC with Indicator; AMP, BAR, COC, MDMA, MET, MTD, OPI, OXY, PCP	E&I Exempt
QODOA-10106EUO-I	10-Panel T-Square with indicator; AMP, BAR, COC, MDMA, MET, MTD, OPI, OXY, PCP, THC	E&I Exempt
QODOA-266-I	6-Panel T-Square with indicator; AMP, BZO, COC, MET, OPI, THC	Forensic Use Only
QODOA-5106	10-Panel T-Square; AMP, BAR, BUP, BZO, COC, MET, OPI, OXY, PCP, THC	Forensic Use Only
QODOA-3126-A*	12-Panel T-Square; AMP, BAR, BUP, BZO, COC, MET, MTD, OPI, OXY, PCP, THC+ALC	Forensic Use Only
QODOA-6126-I	12-Panel T-Square with Indicator; AMP, BAR, BUP, BZO, COC, MDMA, MET, MTD, OPI, OXY, PCP, THC	Forensic Use Only

※ Customizable drug panels & Alcohol test are available.

A* Alcohol testing is included.

Rel.: 2024/03/08