



SAFE life™

Benzodiazepines (BZO) Liquid & Powder Drug Test

- Harm Reduction
- Rapid Results
- Accurate

SAFElife™ Benzodiazepines (BZO) Liquid & Powder Drug Test is an immunochromatographic assay for the qualitative determination of Benzodiazepine in urine, liquid, powder or pill samples.

Benzodiazepines are a class of medications primarily used for their sedative, anxiolytic (anti-anxiety), muscle relaxant, and anticonvulsant properties. They act on the central nervous system (CNS) by enhancing the effect of the neurotransmitter gamma-aminobutyric acid (GABA). In 2017, benzodiazepines and other tranquilizers were the third most commonly used illicit or prescription drug in the U.S. SAFElife™ Benzodiazepine (BZO) Liquid & Powder Drug Test can determine if Benzodiazepine is present in a liquid or powder, providing vital information to help with harm reduction.




Features and Benefits


- Detects Benzodiazepine as well as 17 other analogs, including **Alprazolam, Clonazepam, Diazepam and Lorazepam**.
- No Interference with Fentanyl.
- Clear Interpretation: Results in **1 minute** and are easy to read, with distinct lines indicating the presence of specimen.
- Easy to Use: Each strip features clear **MAX Line** and **MIN Line** indicators.
- Cost-effective option: Test Strip Version and Accessory Kit Version to meet diverse needs
- **All-in-One** Convenience: All the accessories you need in one kit for user-friendly simplicity.

Drug (Identifier)	Benzodiazepine (BZO)
Calibrator	Oxazepam
Cut-off Level	300 ng/mL
Reading Time	1 min
Sample Type	urine, liquid, powder or pill


Detection of Benzodiazepine and Its Analogues				
Oxazepam	Clobazam	Fentanyl	Midazolam	Nordiazepam
Alprazolam	Clonazepam	Flunitrazepam	Nimetazepam	Prazepam
Bromazepam	Diazepam	Lorazepam	Nitrazepam	Temazepam
Chlordiazepoxide	Estazolam	Lormetazepam	Norchlordiazepoxide	Triazolam



Liquid & Powder Drug Test Strip



Liquid & Powder Drug Test Kit Pro



Liquid & Powder Drug Test Kit Premier


Rel.: 2025/06/12




SAFE life™

Test Procedure

Step 1 Add Sample




Step 2 Test



Step 3 Wait for 1min



Step 4 Result



WARNING: THIS TEST CANNOT DETERMINE THE PURITY, COMPOSITION OF URINE, LIQUID, POWDER OR PILL SAMPLES.

Item Code	Description	Packing	Complexity
WSPBZO-101	1-Strip Drug Liquid / Powder Detection Test; BZO Benzodiazepines 300 ng/ml	25T	Forensic Use Only Harm Reduction Use
WSPBZO-101-1PAP	1-Strip Drug Liquid / Powder Detection with Mixing Accessories Pro Version; BZO Benzodiazepines 300 ng/ml	1 Test per Pack	Forensic Use Only Harm Reduction Use