



SAFE life™

T-Swab Multi-Drug Oral Fluid Test



- Tests can be run within 2 minutes
- Screens for 7 common drugs of misuse
- Easy to use
- Result at 5 minutes
- Employment & Insurance Testing

Employment and Insurance (E&I) Testing:

Although validated according to the same rigorous quality system requirements certain Class II devices are exempt from pre-market notification and may be used for Employment and Insurance testing only.

Workplace Drug Testing:

Follow-Up, For Cause, periodic, post-accident, pre-employment, random, return to duty. Not applicable for Federally regulated workplace and DOT testing programs.

T-Swab™ Cutoff Levels

Drug Name	Cut off (ng/ml)
AMP (Amphetamines)	50
COC (Cocaine)	20
THC (Marijuana)	40
MET (Methamphetamine)	50
OPI (Opiates)	40
OXY (Oxycodone)	20
PCP (Phencyclidine)	10

We offers 4 combinations from 7 substances.

Item Code	Test Combination
** SODOA-156NTEUO	AMP-COC-OPI; MET-PCP
** SODOA-256EUO	AMP-COC-OPI; MET-THC
SODOA-166EUO	AMP-COC-OPI; MET-THC-PCP
** SODOA-666EUO	AMP-COC-OPI; MET-THC-OXY

** Ordered upon request

Rel: 2025/08/12