




SAFElife™

T-Swab Multi-Drug Oral Fluid Test

Our compact dual sided all-in-one device delivers 12 drug tests with the added option of alcohol detection.



- Suitable for FOU and Employment & Insurance (E&I) testing scenarios.
- Comprehensive 12-panel drug screening with optional alcohol detection—all from a single, compact device.
- Intuitive design enables accurate testing with minimal training.
- Rapid sample collection and definitive results in 5 minutes.

T-Swab Cutoff Levels

Drug Name	Cutoff Level (ng/mL)	Classification
AMP (Amphetamines)	50	E&I
BAR (Barbiturates)	20	E&I
BUP (Buprenorphine)	5	FUO
BZO (Benzodiazepines)	10	FUO
COC (Cocaine)	20	E&I
MDMA (Methylenedioxymethamphetamine)	50	E&I
MET (Methamphetamine)	50	E&I
MTD (Methadone)	30	E&I
OPI (Opiates)	40	E&I
OXY (Oxycodone)	20	E&I
PCP (Phencyclidine)	10	E&I
THC (Marijuana)	40	E&I
ALC (Alcohol)	>0.02% BAC	FUO

Customizable options to meet your exact needs

Item Code	Test Combination	Classification
SODOA-156NTEUO	AMP-COC-OPI; MET-PCP	E&I
SODOA-256EUO	AMP-COC-OPI; MET-THC	E&I
SODOA-166EUO	AMP-COC-OPI; MET-THC-PCP	E&I
SODOA-666EUO	AMP-COC-OPI; MET-THC-OXY	E&I
SODOA-1106EUO	AMP-COC-OPI; MET-THC-PCP; OXY-MTD; BAR-MDMA	E&I
SODOA-1126	AMP-COC-OPI; MET-THC-PCP; BZO-OXY-MTD; BUP-BAR-MDMA	FUO
SODOA-1136-A	AMP-COC-OPI; MET-THC-PCP; BZO-OXY-MTD; BUP-BAR-MDMA+ALC	FUO

Contact us to place your order.

Rel.: 2025/09/02



+1 (630) 468-2199 || wondfo@wondfousa.com || WondfoUSA.com

Wondfo® || bringing science to life